

CS3 – Bone Cage for Spinal Implants:

Exam Summary

What was the design challenge?

The challenge was to design a spinal fusion cage that supports the spine at first, but gradually transfers load to new bone as it grows — using an adaptive modulus. The aim is to prevent stress shielding and promote long-term fusion.

How does the cage ensure patient safety?

- Ti outer layer on the Apophyseal ring (strong cortical bone) for stability.
- Mg inner layer degrades slowly to reduce stiffness over time (adaptive modulus).
- Chitosan coating on Mg slows degradation, preventing early collapse.
- Slot-in inserts (fast/medium/slow) allow customisation depending on patient bone quality.

How does the design promote bone formation?

- Porous + lattice structure allows blood and nutrient flow.
- Concave pores mimic natural bone pits → osteoblasts love them.
- Nanoscale surface features increase protein adsorption, attract bone cells.
- HA coating (or nanoscale Ti) improves bone bonding.
- rhBMP-2 in centre chamber triggers new bone growth.
- Piezoelectric coatings (like PVDF) create signals that stimulate bone when patient moves.

Why is radiolucency important and how is it achieved?

- Helps doctors see if bone is growing using X-rays/CT.
- Ti causes artefacts, so use Mg or PEEK inner parts (less radiodense).
- Design can be adaptive radiolucent: radiopaque early, radiolucent later as Mg degrades.

Compare different material choices

Material	Pros	Cons
Titanium (Ti)	Strong, good for screw fixation, easy to modify surface	Radiopaque, stress shielding
PEEK	Radiolucent, flexible	Poor cell attachment, hard to shape pores
Mg Alloy	Bioactive, degrades over time, radiolucent	Can degrade too fast (needs coating)
PLA	Biodegradable, proven	Acidic by-products, not ideal for bone growth

What are some added design features?

- Silver coating: Prevents infection, safe in low doses.
- Copper coating: Might help grow blood vessels.
- TPMS lattices: Reduce stiffness but keep strength.
- 3D printed using SLM or EBM, then surface-treated for better bone bonding.

Bonus: A One-Line Summary (for viva or intro)

This spinal cage is a smart, layered implant with a strong Ti outer shell and a resorbable Mg inner core, designed to support spinal fusion safely while stimulating bone growth and allowing doctors to track healing over time.